

Anti-PGBD5 antibody

Catalog Number: 175018

Product name

Anti-PGBD5 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PGBD5

Preparation

This antigen of this antibody was synthetic peptide within human piggybac transposable element-derived protein 5 aa 61-115.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC

Dilutions

WB: 1:10,000-1:20,000

ICC: 1:100-1:200

Validations




Fig1: Western blot analysis on D3 using anti-PGBD5 polyclonal antibody.




Fig2: ICC on HepG2 using anti- PGBD5 polyclonal antibody.