

# ZNF658 antibody



Catalog Number: 117018

## Product name

ZNF658 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

ZNF658 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 71 kDa

## Preparation

This antibody was obtained by immunization of ZNF658 recombinant protein (Accession Number: BC031626). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117018(ZNF658 Antibody) at dilution of 1:1000