

# Cd2ap antibody

Catalog Number: 109018

## Product name

Cd2ap antibody

IP: 1:1000-1:10000

IF: 1:20-1:200

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Cd2ap Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive WB detected in mouse testis tissue, A431 cells, HeLa cells, mouse brain tissue. Positive IP detected in mouse brain tissue, HeLa cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of Cd2ap recombinant protein (Accession Number: NM\_009847). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109018(Cd2ap antibody) at dilution of 1:200



IP Result of anti-Cd2ap (IP:Catalog No:109018, 4ug; Detection:Catalog No:109018 1:2000) with mouse brain tissue lysate 2560ug.



Immunofluorescent analysis of A431 cells using Catalog No:109018(Cd2ap Antibody) at dilution of

# Cd2ap antibody



Catalog Number: 109018

---

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)