

Human MTHFS (His tag) recombinant protein



Catalog Number: 519018

General Information

Gene Name Synonym

5-formyltetrahydrofolate cyclo-ligase,5,10-methenyl-tetrahydrofolate synthetase,MTHFS,Methenyl-THF synthetase

Protein Construction

A DNA sequence encoding the Human MTHFS Met1-Ala203 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris,200mM NaCl,1mM DTT,50% Glycerol,pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human MTHFS (His tag) recombinant protein