

ORAI3 antibody



Catalog Number: 113418

Product name

ORAI3 antibody

Immunogen

[Human ORAI3 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ORAI3 Rabbit Polyclonal antibody. Positive IHC detected in human endometrial cancer tissue. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: ~30 kDa

Preparation

This antibody was obtained by immunization of ORAI3 recombinant protein (Accession Number: NM_152288). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113418(ORAI3 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using Catalog No:113418(ORAI3 Antibody) at dilution of 1:50 (under 10x lens)



ORAI3 antibody



Catalog Number: 113418

Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using

Catalog No:113418(ORAI3 Antibody) at dilution of 1:50 (under 40x lens)