

FBRS Polyclonal Antibody



Catalog Number: 164418

Product name

FBRS Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to FBRS

Preparation

Antigen: Recombinant protein of human FBRS

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:200 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:100

Validations



Western blot - FBRS Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using FBRS antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.