

Anti-CD1E antibody



Catalog Number: 176418

Product name

WB, IHC-P

Anti-CD1E antibody

Dilutions

Specificity

WB: 1:500-1:2000

Human

IHC-P: 1:100-300

Antibody description

Rabbit polyclonal antibody to CD1E

Validations

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human cd1e. aa range:321-370

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications




Fig1:: Western Blot analysis of K562 cells using CD1E Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:100