

TRPV5 antibody



Catalog Number: 116418

Product name

TRPV5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRPV5 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells.

Preparation

This antibody was obtained by immunization of TRPV5 recombinant protein (Accession Number: BC034740). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HEK-293 cells using Catalog No:116418 (TRPV5 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG