

Anti-NFkB p100/p52 Rabbit antibody



Catalog Number: 167418

Product name

Anti-NFkB p100/p52 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NFkB p100/p52 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human NFkB p100 / p52

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

IP: 1/50

Validations



Western blot analysis of NFkB p100 / p52 expression in (1) Jurkat cell lysate; (2) Mouse heart lysate.