

# c-MYC antibody

Catalog Number: 109418

## Product name

c-MYC antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

## Specificity


Human, Mouse; other species not tested.

IHC: 1:20-1:200

## Antibody description

c-MYC Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human breast cancer tissue, human prostate cancer tissue. Positive IP detected in HL-60 cells. Positive WB detected in Raji cells, HL-60 cells, mouse embryo tissue. Observed molecular weight by Western-blot: 62-65 kDa, 50 kDa

## Validations




## Preparation

This antibody was obtained by immunization of c-MYC recombinant protein (Accession Number: NM\_002467). Purification method: Antigen affinity purified.

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109418(MYC antibody) at dilution of 1:1000

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:109418(MYC antibody) at dilution of 1:50 (under 40x lens)

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP, IHC



## Dilutions

Recommended Dilution:

# c-MYC antibody



Catalog Number: 109418

---

IP Result of anti-MYC (IP:Catalog No:109418, 3ug; Detection:Catalog No:109418 1:500) with HL-60 cells lysate 1200ug.