

GAD65 antibody

Catalog Number: 110818

Product name

GAD65 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GAD65 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000818). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IF: 1:10-1:100


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110818(GAD2 antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:110818(GAD2 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



IP Result of anti-GAD2 (IP:Catalog No:110818, 4ug; Detection:Catalog No:110818 1:1000) with mouse brain tissue lysate 3440ug.