

# GPR45 Polyclonal Antibody



Catalog Number: 162818

## Product name

GPR45 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to GPR45

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human GPR45 (NP\_009158.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - GPR45 Polyclonal Antibody

Western blot analysis of extracts of SH-SY5Y cells, using GPR45 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.