

# Anti-EPHA1 antibody



Catalog Number: 175818

## Product name

Anti-EPHA1 antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Rat, Mouse

## Applications

WB, IHC-P

## Antibody description

Rabbit polyclonal antibody to EPHA1

## Dilutions

WB: 1:1000-1:2000

## Preparation

This antigen of this antibody was synthetic peptide of epha1

IHC-P: 1:100-200

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal




Fig1:: Western blot analysis of 1) Mouse Brain Tissue, 2) Hela, 3) Rat Brain Tissue with EphA1 Rabbit pAb diluted at 1:2,000.