

CALCOCO2 antibody



Catalog Number: 108818

Product name

CALCOCO2 antibody

Immunogen

[Human CALCOCO2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CALCOCO2 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue, human skeletal muscle tissue. Positive WB detected in HeLa cells, human placenta tissue, Jurkat cells, mouse skeletal muscle tissue, rat skeletal muscle tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of CALCOCO2 recombinant protein (Accession Number: NM_005831). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108818(CALCOCO2 antibody) at dilution of




1:1000

IP Result of anti-CALCOCO2 (IP:Catalog No:108818, 3ug; Detection:Catalog No:108818 1:1000) with HeLa cells lysate 2150ug.


CALCOCO2 antibody



Catalog Number: 108818



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108818(CALCocO2 antibody) at dilution of



1:100 (under 10x lens)

rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108818(CALCocO2 Antibody) at dilution of 1:1000