

General Information

Protein Construction

A DNA sequence encoding the ferret IFNG (A4PIZ9) (Met1-Lys166) was expressed with a polyhistidine tag at the C-terminus.

Organism

Ferret

Expression Host

Human Cells

QC Testing

Purity

(43.0+47.1+9.3) % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 24

Molecular Mass

The recombinant Ferret IFNG consists of 154 amino acids and has a calculated molecular mass of 18.4 kDa. The apparent molecular mass of the recombinant protein is approximately 28,23 and

19 KDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Ferret IFNG / Interferon Gamma Protein (His Tag)
SDS-PAGE