

Netrin G1 antibody



Catalog Number: 113118

Product name

Netrin G1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human,Mouse,Rat; other species not tested.

Applications

ELISA, WB

Antibody description

Netrin G1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, A549 cells, human brain tissue, human testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 60 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Netrin G1 recombinant protein (Accession Number: NM_001372168). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:113118(NTNG1 antibody) at dilution of 1:500