

PPCS antibody



Catalog Number: 114118

Product name

PPCS antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPCS Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human spleen tissue, human testis tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of PPCS recombinant protein (Accession Number: NM_024664). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114118(PPCS antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114118(PPCS antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114118(PPCS antibody) at dilution of 1:100 (under 40x lens)