

Anti-NEK5 antibody



Catalog Number: 176118

Product name

Anti-NEK5 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to NEK5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human nek5. aa range:160-209

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using Nek5 Polyclonal Antibody