

# TIA1 antibody

Catalog Number: 116118

## Product name

TIA1 antibody

## Immunogen

[Human TIA1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TIA1 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, mouse spleen tissue, mouse thymus tissue, Raji cells. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of TIA1 recombinant protein (Accession Number: NM\_001351518). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116118(TIA1 antibody) at dilution of 1:200



IP Result of anti-TIA1 (IP:Catalog No:116118, 4ug; Detection:Catalog No:116118 1:500) with Jurkat cells lysate 3200ug.