

Anti-eIF4A1 Rabbit antibody



Catalog Number: 167118

Product name

Anti-eIF4A1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

eIF4A1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human eIF4A1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of eIF4A1 in K562,C6,3T3,Hela cell lysates using eIF4A1 Rabbit pAb(1:1000 diluted).Predicted band size:46kDa.Observed band size:46kDa.