

# Brentuximab (Synonyms: MOAB, CHIMERIC, SGN-30, TO CD30 ANTIGEN)



Catalog Number: 177118

## Product name

Brentuximab (Synonyms: MOAB, CHIMERIC, SGN-30, TO CD30 ANTIGEN)

## Specificity

Human TNFRSF8 / CD30

## Antibody description

Anti-TNFRSF8 / CD30 Antibody (Brentuximab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-TNFRSF8 / CD30 Antibody (brentuximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFRSF8 / CD30 Antibody (brentuximab) is 98.99%, determined by SEC-



## HPLC.

Bioactivity: ELISA

Immobilized human CD30 His at 2 ug/mL can bind Anti-TNFRSF8 / CD30 Antibody (brentuximab), EC50=0.006007 ug/mL.