

Anti-CBX1 Mouse antibody

Catalog Number: 169118

Product name

Anti-CBX1 Mouse antibody

Specificity

Human, Monkey

Antibody description

CBX1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human PAX3 (AA: 142-203) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB, ICC/IF, IHC-P, FC

Validations




Flow cytometric analysis of COS7 cells using CBX1 mouse mAb (green) and negative control



(purple). 17-

Western blot analysis using CBX1 mouse mAb against HeLa (1), COS7 (2), NIH/3T3 (3), A431 (4), and C6 (5) cell lysate.



Immunofluorescence analysis of HeLa cells using CBX1 mouse mAb (green). Blue