

HOMER2 antibody



Catalog Number: 111518

Product name

HOMER2 antibody

Specificity

Human; other species not tested.

Antibody description

HOMER2 Rabbit Polyclonal antibody. Positive WB detected in U2OS cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of HOMER2 recombinant protein (Accession Number: NM_004839). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



WB result of anti-HOMER2 (Catalog No:111518 1:2000) with transfected U2OS cell by Dr. Yupeng Chen.