

S100A11 antibody



Catalog Number: 107518

Product name

S100A11 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

S100A11 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, A431 cells, human lung tissue. Positive IHC detected in human skin. Observed molecular weight by Western-blot: 12 kDa

Preparation

This antibody was obtained by immunization of S100A11 recombinant protein (Accession Number: NM_005620). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107518(S100A11 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human skin using Catalog No:107518(S100A11 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human skin using Catalog No:107518(S100A11 antibody) at dilution of 1:50 (under 40x lens)