

CPNE6 antibody



Catalog Number: 109518

Product name

CPNE6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CPNE6 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of CPNE6 recombinant protein (Accession Number: XM_005268216). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109518(CPNE6 antibody) at dilution of 1:500