

kazrin antibody

Catalog Number: 111918

Product name

kazrin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

kazrin Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in human heart tissue, mouse brain tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of kazrin recombinant protein (Accession Number: NM_001018000). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111918(RP1-21O18.1 antibody) at dilution of 1:1200



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111918(RP1-21O18.1 antibody) at dilution of 1:100 (under 10x lens)