

# ZBTB46 antibody



Catalog Number: 116918

## Product name

ZBTB46 antibody

## Specificity

Human, Rat, Mouse; other species not tested.

## Antibody description

ZBTB46 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of ZBTB46 recombinant protein (Accession Number: NM\_001369741). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116918(ZBTB46 Antibody) at dilution of 1:600