

AGER antibody

Catalog Number: 107918

Product name

AGER antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AGER Rabbit Polyclonal antibody. Positive IHC detected in human renal cell carcinoma tissue, human lung cancer tissue. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of AGER recombinant protein (Accession Number: NM_001136). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107918(AGER Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human renal cell carcinoma tissue slide using Catalog No:107918(AGER Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human renal cell carcinoma tissue slide using

AGER antibody



Catalog Number: 107918

Catalog No:107918(AGER Antibody) at dilution of 1:200 (under 40x lens).