

Anti-ATP6V1A Rabbit antibody



Catalog Number: 168918

Product name

Anti-ATP6V1A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ATP6V1A Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human ATP6V1A

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

1:1000


1:200

1:50

1:20

1:20

Validations



Western blot detection of ATP6V1A in K562, Rat Brain, 3T3 cell lysates using ATP6V1A Rabbit mAb(1:1000 diluted). Predicted band size: 68kDa. Observed band size: 68kDa.