

Growth hormone antibody



Catalog Number: 111218

Product name

Growth hormone antibody

Specificity

Human; other species not tested.

Antibody description

Growth hormone Rabbit Polyclonal antibody.
Positive WB detected in human placenta tissue.
Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of Growth hormone recombinant protein (Accession Number: NM_000515). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:111218(GH1 antibody) at dilution of 1:1000