

LHX5 antibody



Catalog Number: 112218

Product name

LHX5 antibody

Specificity

Human; other species not tested.

Antibody description

LHX5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, human brain tissue. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of LHX5 recombinant protein (Accession Number: NM_022363). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112218(LHX5 antibody) at dilution of 1:300