

# Anti-PTPMT1 antibody



Catalog Number: 106218

## Product name

Anti-PTPMT1 antibody

## Immunogen

[Human PTPMT1 \(His Tag\) recombinant protein](#)

## Specificity

Human PTPMT1

## Antibody description

Rabbit polyclonal to PTPMT1

## Preparation

Produced in rabbits immunized with purified, recombinant Human PTPMT1 (rh PTPMT1; Q8WUK0-1; Lys28-Thr201). PTPMT1 specific IgG was purified by Human PTPMT1 affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 2-5  $\mu$ g/ml


IHC-P: 0.1-2  $\mu$ g/ml

## Validations




Items	Lanes	A
Sample (whole cell lysate)		293T
Sample Volume ( $\mu$ g/lane)		60
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		2.5 $\mu$ g/ml
Secondary Antibody		Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution

PTPMT1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



PTPMT1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry




# Anti-PTPMT1 antibody

Catalog Number: 106218

---

PTPMT1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry



PTPMT1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry