

Anti-DDIT3 Rabbit antibody



Catalog Number: 167218

Product name

Anti-DDIT3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

DDIT3 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human DDIT3

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/500-1/2000

IHC: 1/50-1/200

FC: 1/50

Validations



Western blot analysis of DDIT3 expression in (1) HeLa cell lysate; (2) NIH/3T3 cell lysate; (3) C6 cell lysate.