

# Anti-S100A4 (3D10) Mouse antibody



Catalog Number: 168218

## Product name

Anti-S100A4 (3D10) Mouse antibody

## Specificity

Human

## Antibody description

S100A4 (3D10) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant human S100A4 protein.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, IP


## Dilutions

WB: 1/500

## Validations



Immunoprecipitation analysis of HeLa cell lysates



using S100A4 mouse mAb.

Western blot detection of S100A4 in HeLa and A549 cell lysates using S100A4 mouse mAb(dilution 1:500).Predicted band size:12kDa.Observed band size:12kDa.