

# SSX2IP antibody



Catalog Number: 115618

---

## Product name

SSX2IP antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SSX2IP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of SSX2IP recombinant protein (Accession Number: XM\_047444216). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115618(SSX2IP antibody) at dilution of 1:500