

# UNC5C-Specific antibody



Catalog Number: 116618

---

## Product name

UNC5C-Specific antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

UNC5C-Specific Rabbit Polyclonal antibody.  
Positive WB detected in mouse brain tissue.  
Observed molecular weight by Western-blot:  
50kd,103kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003728).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116618(UNC5C-Specific antibody) at dilution of 1:1000