

Anti-Glucocorticoid Receptor Rabbit antibody



Catalog Number: 167618

Product name

Anti-Glucocorticoid Receptor Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Glucocorticoid Receptor Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Glucocorticoid Receptor

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/2000-1/10000

FC: 1/20-1/100

Validations



Western blot detection of Glucocorticoid Receptor in CHO-K1, C6, Rat Brain, Ramos cell lysates using Glucocorticoid Receptor Rabbit mAb(1:1000 diluted). Predicted band size: 86kDa. Observed band size: 86kDa.