

C19orf20 antibody



Catalog Number: 108618

Product name

C19orf20 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C19orf20 Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue, mouse heart tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 29KD

Preparation

This antibody was obtained by immunization of C19orf20 recombinant protein (Accession Number: BC009520). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:108618(C19orf20 Antibody) at dilution of 1:600