

Anti-MEP1A antibody



Catalog Number: 102318

Product name

Anti-MEP1A antibody

Immunogen

[Mouse MEP1A \(His Tag\) recombinant protein](#)

Specificity

Mouse MEP1A

Antibody description

Rabbit polyclonal to MEP1A

Preparation

Produced in rabbits immunized with purified, recombinant Mouse MEP1A (rM MEP1A; NP_032611.2; Met 1-Arg 615). Total IgG was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 10-20 µg/ml

ELISA: 0.5-1.0 µg/mL

This antibody can be used at 0.5-1.0 µg/mL with the appropriate secondary reagents to detect Mouse MEP1A. The detection limit for Mouse MEP1A is 0.00245 ng/well.

Validations



Lanes	A
Items	
Sample (Membrane, lysate)	A431
Sample Volume (µg/lane)	20
Gel	12% SDS-PAGE reducing gel
Recommended Concentration	10-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG [H+], at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 84 kDa Observed band size : 70 kDa

MEP1A / Endopeptidase-2 Antibody, Rabbit PAb, Western blot