

# SLC24A4 antibody



Catalog Number: 115318

## Product name

SLC24A4 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

SLC24A4 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 62-90kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_153646). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115318(SLC24A4 antibody) at dilution of 1:300