

TRIM46 antibody



Catalog Number: 116318

Product name

TRIM46 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TRIM46 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, HeLa cells, mouse testis tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 55kd,62kd

Preparation

This antibody was obtained by immunization of TRIM46 recombinant protein (Accession Number: BC069416). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116318(TRIM46 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using TRIM46 antibody Catalog No:116318 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).