

Anti-CDA Rabbit antibody



Catalog Number: 167318

Product name

Anti-CDA Rabbit antibody

Specificity

Human

Antibody description

CDA Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human CDA

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/1000

IHC: 1/500

FC: 1/20

Validations



Western blot detection of CDA in HeLa cell lysates using CDA Rabbit pAb(1:1000 diluted).Predicted band size:16kDa.Observed band size:16kDa.