

# CIP29 antibody



Catalog Number: 109318

---

## Product name

CIP29 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CIP29 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HEK-293 cells, HeLa cells, HepG2 cells, Raji cells, rat liver tissue, rat pancreas tissue. Observed molecular weight by Western-blot: 29-32 kDa

## Preparation

This antibody was obtained by immunization of CIP29 recombinant protein (Accession Number: NM\_033082). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109318(CIP29 antibody) at dilution of 1:1500