

FPGT antibody

Catalog Number: 110718

Product name

FPGT antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FPGT Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue. Positive IHC detected in human colon tissue, human kidney tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of FPGT recombinant protein (Accession Number: NM_003838). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:110718(FPGT Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110718(FPGT Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110718(FPGT Antibody) at dilution of 1:50

FPGT antibody



Catalog Number: 110718

(under 40x lens)