

Anti-CKM antibody



Catalog Number: 142718

Product name

Anti-CKM antibody

Immunogen

[Human CKM / CK-MM \(His Tag\) recombinant protein](#)

Specificity

Human CKM / CK-MM

Antibody description

Rabbit Polyclonal to Human CKM

Preparation

Produced in rabbits immunized with purified, recombinant Human CKM / CK-MM . CKM / CK-MM specific IgG was purified by Human CKM / CK-MM affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human CKM / CK-MM.

IHC-P: 0.1-2 µg/mL

Validations



Human CKM / CK-MM
Immunohistochemistry(IHC) 14730