

ACADL antibody



Catalog Number: 107718

Product name

ACADL antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

ACADL Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells, HeLa cells, HepG2 cells, NIH/3T3 cells. Positive IHC detected in human kidney tissue, human liver tissue. Observed molecular weight by Western-blot: 45-48 kDa

Preparation

This antibody was obtained by immunization of ACADL recombinant protein (Accession Number: NM_001608). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107718(ACADL antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107718(ACADL antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107718(ACADL antibody) at dilution of 1:50 (under 40x lens)