

C4 (gamma chain) antibody



Catalog Number: 107028

Product name

C4 (gamma chain) antibody

Specificity

Human; other species not tested.

Antibody description

C4 (gamma chain) Mouse Monoclonal antibody.
Positive WB detected in human plasma tissue.
Observed molecular weight by Western-blot: 31 kDa

Preparation

This antibody was obtained by immunization of C4 (gamma chain) recombinant protein (Accession Number: BC012372). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107028(C4A Antibody) at dilution of 1:8000