

Anti-ZNF364/RFN115 Rabbit antibody



Catalog Number: 169028

Product name

Anti-ZNF364/RFN115 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

ZNF364/RFN115 Rabbit monoclonal antibody

Dilutions

1:1000


Preparation

Antigen: A synthetic peptide of human ZNF364

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of ZNF364/RFN115 in Jurkat, HeLa cell lysates using ZNF364/RFN115 Rabbit mAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.

Clonality

Monoclonal