

OSBP antibody



Catalog Number: 113428

Product name

OSBP antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

OSBP Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, mouse kidney tissue, rat brain tissue, rat kidney tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 89kd

Preparation

This antibody was obtained by immunization of OSBP recombinant protein (Accession Number: NM_002556). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions

Recommended Dilution:

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113428(OSBP



antibody) at dilution of 1:800

IP result of anti-OSBP(Catalog No:113428 for IP and Detection).



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog

OSBP antibody



Catalog Number: 113428

No:113428(OSBP antibody) at dilution of 1:50 (under 10x lens)