

# OSBP antibody

Catalog Number: 113428

## Product name

OSBP antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

OSBP Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, mouse kidney tissue, rat brain tissue, rat kidney tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 89kd

## Validations




## Preparation

This antibody was obtained by immunization of OSBP recombinant protein (Accession Number: NM\_002556). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113428(OSBP

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

antibody) at dilution of 1:800

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

IP result of anti-OSBP(Catalog No:113428 for IP and Detection).



## Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog

# OSBP antibody



Catalog Number: 113428

---

No:113428(OSBP antibody) at dilution of 1:50 (under 10x lens)